

158
60

INDIAN INSTITUTE
OF
TECHNOLOGY
PHYSICS DEPT.

NO. 190100079 CLASS

BATCH P19/1

NAME NABH AGRAWAL

EXPT. No. 1

LABORATORY

DATE 05/08/2019

AIM

To observe the intensity pattern generated by diffraction of a laser beam when passed through a double slit and to measure the size of diffracting element.

APPARATUS

- i) Optical Rail.
- ii) Diode Laser (wavelength 650 nm).
- iii) Kinematic Laser Mount.
- iv) Power Supply for LASER.
- v) Cell Mount with Diffraction Cell.
- vi) Pin hole detector.
- vii) Output Measurement Unit.
- viii) Linear Translational Stage
- ix) Slit Box.

THEORY & FORMULAE USED

Light rays, when partially obstructed tend to show certain wave like properties that include diffraction and interference (in case of overlapping of waves). This bending of light was explained by Huygen's Principle which treated every point on wave front

as another source of Light.

1) For a diffraction by Single Slit, intensity is given by

$$I(\theta) = I_0 \left(\frac{\sin \beta}{\beta} \right)^2 \quad \text{where } \beta = \left(\frac{\pi d}{\lambda} \right) \sin \theta.$$

Thus, for a maxima, $\beta = 0 \& \pm (2m+1)\frac{\pi}{2}$; $m=1, 2, \dots$
 $\theta = 0 \& \pm \sin^{-1} \left[(2m+1)\frac{1}{2d} \right]$.

for a minima, $\beta_m = \pm m\pi$; $m=1, 2, \dots$
 $\theta_m = \pm \sin^{-1} \left(\frac{m\pi}{d} \right)$.

∴ From the positions of minima, we can find out the slit width using $[d = m\lambda / \sin \theta_m]$, where θ_m could be experimentally calculated using $\boxed{\tan \theta_m = \frac{y_m}{D}}$.

y_m = y-position of Minima.

D = distance between Slit & screen.

d = Slit width.

λ = wavelength used.

2) In Case of Interference & Diffraction by Double Slit:

$$I(\theta) = I_0 \left(\frac{\cos^2 \left(\frac{\pi d \sin \theta}{\lambda} \right)}{1 + \left(\frac{\pi d \sin \theta}{\lambda} \right)^2} \right) \quad \text{where } \beta = \left(\frac{\pi d}{\lambda} \right) \sin \theta.$$

a = Slit width

d = Slit separation.

∴ for a maxima $d \sin \theta = m\lambda$, $m \in \mathbb{Z}$.

for a minima $d \sin \theta = (m + \frac{1}{2})\lambda$, $m \in \mathbb{Z}$.

Here $d \gg a$, and from positions of maxima, we can easily find the slit separation: $\boxed{d = \frac{m\lambda}{\sin \theta_m}}$ where $\boxed{\theta_m = \tan^{-1} \left(\frac{y_m}{D} \right)}$

INDIAN INSTITUTE
OF
TECHNOLOGY

PHYSICS DEPT.

NO. (90100079) CLASS

BATCH P14/1

NAME NABH AGRAWAL

EXPT. No. 1

DATE 05/08/2019

LABORATORY

3) Diffraction of Light by Pin hole.

$$I(\theta) = I_0 \left[\frac{2 J_1(k a \sin \theta)}{k a \sin \theta} \right]^2 \quad \text{where } k = 2\pi/d \quad \& \quad J_1 = \text{Bessel Function of first kind of order one.}$$

Bright rings appear around the central maxima (called Airy Disc)


Minima occurs at zeros of J_1 at $k a \sin \theta = 3, 83, 7.0156, \dots$
And for first minima, $\sin \theta = 1.22 d/a$

$$\text{Least Count} = (\text{Smallest Division on Main Scale}) / (\text{No. of Division on Circular Scale.})$$

$$= 0.05 \text{ cm} / 50 = 0.001 \text{ cm.}$$

$$\text{Micrometer Reading} = \text{M.S.R.} + (\text{No. of Div. on Circular Scale}) \times \text{L.C.}$$

DIAGRAM



"DOUBLE SLIT SYSTEM"

OBSERVATIONS

Difraction Pattern by Double Slit.

Main Scale Reading	Circular Scale Reading	Total Reading (mm)	Intensity (μA)
18.0	0	18.0	09.0
16.5	0	16.5	04.4
15.0	0	15.0	09.1
15.5	0	15.5	05.1
15.0	0	15.0	03.8
14.5	0	14.5	03.5
14.0	0	14.0	03.7
13.5	0	13.5	08.5
13.0	32	13.32	09.7
13.0	0	13.0	05.8
12.5	20	12.70	04.0
12.0	0	12.5	03.5
12.0	0	12.0	05.0
11.5	0	11.5	20.1
11.0	40	11.40	39.6
11.0	30	11.30	62.6
11.0	25	11.25	75.1
11.0	20	11.20	85.9
11.0	5	11.05	99.2
11.0	0	11.0	96.3
10.5	20	10.70	29.2
10.5	0	10.5	52.2
10.0	30	10.30	158.5

INDIAN INSTITUTE
OF
TECHNOLOGY

PHYSICS DEPT.

NO. 190100079 CLASS

BATCH P14/1

NAME NABH AGRAWAL

EXPT. No. 1

LABORATORY

DATE

05/08/2019

Main Scale Reading	Circular Scale Reading	Total Reading (mm)	Intensity (uA)
10.0	0	10.0	172.0
9.5	5	9.55	30.0
9.5	0	9.5	37.2
9.0	15	9.15	82.5
9.0	0	9.0	66.3
8.0	0	8.0	03.0
7.5	0	7.5	03.2
7.0	0	7.0	11.3
6.5	40	6.90	11.8
6.5	0	6.50	04.2
6.0	10	6.10	02.1
5.5	0	5.5	02.3
5.0	15	5.15	02.4
5.0	5	5.05	03.2
5.0	0	5.0	03.6
4.5	25	4.75	05.0
4.5	0	4.5	03.7
4.0	0	4.0	02.3

(a) ✓

$$D = 1100 \text{ mm} - 435 \text{ mm} = 665 \text{ mm}$$

5.8 min
05/08/2019

PRECAUTIONS

- 1) Do not touch the inner surface of diffraction cell.
- 2) Avoid Backlash error while moving the micrometer scale.
- 3) Take care while handling all the equipments.
- 4) Avoid looking directly into the laser beam.

CALCULATIONS

$$y_1 \text{ Position of Central Maxima} = 10.00 \text{ mm}$$

$$\text{Positions of } 1^{\text{st}} \text{ Maxima} = 11.05 \text{ mm} \text{ & } 9.15 \text{ mm}$$

$$\text{Positions of } 2^{\text{nd}} \text{ Maxima} = 13.32 \text{ mm} \text{ & } 6.90 \text{ mm}$$

$$\text{Positions of } 3^{\text{rd}} \text{ Maxima} = 15.50 \text{ mm} \text{ & } 4.75 \text{ mm}$$

$$y_1 = 11.05 - 10.00 = 1.05 \text{ mm}; y'_1 = 10.00 - 9.15 = 0.85 \text{ mm}$$

$$y_2 = 13.32 - 10.00 = 3.32 \text{ mm}; y'_2 = 10.00 - 6.90 = 3.10 \text{ mm}$$

$$y_3 = 15.50 - 10.00 = 5.50 \text{ mm}; y'_3 = 10.00 - 4.75 = 5.25 \text{ mm}$$

$$\text{Also, } D = 665 \text{ mm.}$$

Now, Using $\theta_m = \tan^{-1}\left(\frac{y_m}{D}\right)$. we have,

$$\theta_1 = 0.0904; \quad ; \quad \theta'_1 = 0.0732.$$

$$\theta_2 = 0.2860; \quad ; \quad \theta'_2 = 0.2670.$$

$$\theta_3 = 0.4738; \quad ; \quad \theta'_3 = 0.4523.$$

INDIAN INSTITUTE
OF
TECHNOLOGY
PHYSICS DEPT.

NO. 190100079 CLASS

BATCH P14/1

NAME NABH AGRAWAL

EXPT. No. 1

LABORATORY

DATE 05/08/2019

Now using $d = \frac{m d}{\sin \Theta_m}$ we have,

$$d_1 = 411.972 \mu\text{m} ; \quad d_4 = 508.774 \mu\text{m.} \quad \checkmark$$

$$d_2 = 260.436 \mu\text{m.} ; \quad d_5 = 278.969 \mu\text{m.} \quad \checkmark$$

$$d_3 = 235.812 \mu\text{m.} ; \quad d_6 = 247.021 \mu\text{m.} \quad \checkmark$$

$$\therefore d_{avg} = \frac{\sum d_i}{6} = \frac{411.972 + 508.774 + 260.436 + 278.969 + 235.812 + 247.021}{6} \quad \checkmark$$

$$= 323.830 \mu\text{m.} \quad \checkmark$$

$$\begin{aligned} \delta d &= \sqrt{\frac{\sum (d_i - d_{avg})^2}{6}} = \sqrt{\frac{7769.0121 + 34204.2831 + 4018.7982 + 2012.5093 + 7747.1683 + 5899.6224}{6}} \\ &= \sqrt{10275.2324}, \\ &= 101.3668 \mu\text{m} \quad \checkmark \end{aligned}$$

10

RESULT

$$\text{Slit Width} = d_{\text{avg}} \pm \delta d$$

$$= 323.830 \mu\text{m} \pm 101.3668 \mu\text{m}$$

Scanned by CamScanner

Diffraction Pattern by Double Slit

SCALE :-
 x: 1 Smallest Box = 0.1 mm.
 y: 1 Smallest Box = 1 μA.

